

## FOCUSED TRI-POLAR ALTERNATING DEPTH RF



# YOUR CLIENTS' SKIN

- Skin Tightening
- Body Shaping, Slimming
- · Collagen Remodeling
- Cellulite Reduction





## FOCUSED TRI-POLAR ALTERNATING DEPTH RF

Non-invasive, non-surgical facial and body treatments are among the largest and fastest growing areas in aesthetic medicine today. The **TRIO** XT, is the newest radio frequency (**RF**) system from **Sybaritic**. The **TRIO** XT offers both auto-focus bi-polar and the new proprietary advanced technology feature called "**Focused Alternating Depth (FAD**\*\*) **TRI-Polar RF**".

Focused Alternating Depth technology (patent pending) is incorporated in the **TRIO** XT large body **TRI-Polar** handpiece for greater **RF** energy coverage and more efficient penetration. **FAD** reaches both the dermal-epidermal junction and the dermal-subdermal junction, improving body contouring and cellulite reduction results.

FAD TRI-Polar technology delivers a continuous RF signal that focuses into two different depths in the treatment area during the same cycle. This dual depth penetration is a significant improvement because it sends two alternating signals (to two levels of the treatment area) instead of only one signal to a single depth level. FAD provides denser energy into the target area and vertically penetrates to both medium and deeper depth levels for more thorough vertical heating and faster denaturation. The results show immediate skin tightening, reduction of circumferential diameter and diminished fine lines and wrinkles. FAD also helps you achieve better tone and texture as a result of more effective cellulite reduction. The return on investment is extremely fast, which makes the TRIO XT the best value for your practice.

#### **TRIO XT Features**

- Focused TRI-Polar Alternating Depth (FAD) Technology
- Ten (10) visually displayed treatment programs
- Large 25.4 cm touch screen adjustable LCD display
- Bi-Polar and **FAD TRI-Polar** radio frequency ensure precise electro-thermolysis
- Highly effective Clinical Instruction Guide

#### Three (3) Versatile Treatment Handpieces:

- 1. TRI-Polar FAD Body Handpiece
- 2. Bi-Polar Variable Contact Cooled Facial Handpiece
- 3. Bi-Polar Fine Tuning Mini-Facial Probe Handpiece



# Tightens, Tones & Redefines Radiant Beautiful Skin



## TECHNICAL SPECIFICATIONS

**Technology:** Focused Alternating Depth (Patent Pending)

Programs: Ten (10) Pre-set

**Color LCD Touch screen:** Large 25.4 cm / 10" **RF Carrier Frequency (CF):** Low Frequency 5 MHz

**RF Modulating Frequencies (MF):** Proprietary Ultra-Low Frequency **RF Proprietary Sequence Modulation:** 1 Hz, 2 Hz, 3 Hz, 5 Hz, 10 Hz

Energy Adjustment: 0 - 100%

Contact Cooling: Three (3) Temperatures Selections, -5° C, 2.5° C 10° C

RF Power: Proprietary, Adjustable, 50 W

Facial Handpiece: 2.5 cm (with contact cooling and Bi-Polar)

Mini-Probe Handpiece: 2 cm (Bi-Polar)

**Body Handpiece:** 9.6 cm (Alternating Depth TRI-Polar)

**Dimensions:** 15" L x 21" W x 51" H; .38 m L x .53 m W x 1.29 m H

Device Weight: 132 lb. / 60 kg.

Electrical Requirement: 110 V or 220 V, 50-60 Hz

Extended Warranty Available EU/CE Certifications Granted





SYBARITIC INC., Global Aesthetic Solutions
9220 James Avenue South Minneapolis, Minnesota 55431 USA · www.sybaritic.com